

bs-21228R

[Primary Antibody]

SLC25A30 Rabbit pAb

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 253512</p> <p>Target: SLC25A30</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SLC25A30: 1-100/291.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Although the outer mitochondrial membrane is permeable to many small metabolites, transport of solutes across the inner mitochondrial membrane is achieved by members of the mitochondrial carrier protein family, such as SLC25A30 (Haguenauer et al., 2005 [PubMed 15809292]).[supplied by OMIM, Mar 2008]</p>	<p>Isotype: IgG</p> <p>SWISS: Q5SVS4</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat)</p> <p>Predicted MW.: 32 kDa</p> <p>Subcellular Location: Cytoplasm</p>
---	--	---