

bs-21232R

[Primary Antibody]

SLC25A34 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 284723</p> <p>Target: SLC25A34</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SLC25A34: 21-120/304.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: SLC25A34 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).[supplied by OMIM, Mar 2008]</p>	<p>Isotype: IgG</p> <p>SWISS: Q6PIV7</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat)</p> <p>Predicted MW.: 32 kDa</p> <p>Subcellular Location: Cytoplasm</p>
--	--	---