
SLC25A43 Rabbit pAb

Catalog Number: bs-21239R

Target Protein: SLC25A43

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Sheep, Cow, Chicken, Dog)

Predicted MW: 38 kDa

Entrez Gene: 203427

Swiss Prot: Q8WUT9

Source: KLH conjugated synthetic peptide derived from human SLC25A43: 251-350/341.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a member of the mitochondrial carrier family of proteins.[provided by RefSeq, Dec 2009]