

bs-21239R

[Primary Antibody]

SLC25A43 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|----------------------|--|
| Host: Rabbit | Isotype: IgG | Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Chicken, Dog) Predicted MW.: 38 kDa Subcellular Location: Cytoplasm |
| Clonality: Polyclonal | | |
| GeneID: 203427 | SWISS: Q8WUT9 | |
| Target: SLC25A43 | | |
| Immunogen: KLH conjugated synthetic peptide derived from human SLC25A43: 251-350/341. | | |
| Purification: affinity purified by Protein A | | |
| Concentration: 1mg/ml | | |
| Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. | | |
| Background: This gene encodes a member of the mitochondrial carrier family of proteins.[provided by RefSeq, Dec 2009] | | |