

bs-21239R

[Primary Antibody]

SLC25A43 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 203427</p> <p>Target: SLC25A43</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SLC25A43: 251-350/341.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene encodes a member of the mitochondrial carrier family of proteins.[provided by RefSeq, Dec 2009]</p>	<p>Isotype: IgG</p> <p>SWISS: Q8WUT9</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Chicken, Dog)</p> <p>Predicted MW.: 38 kDa</p> <p>Subcellular Location: Cytoplasm</p>
--	--	---