bs-9182R

[Primary Antibody]

Bioss ANTIBODIES

VGLL1 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 51442 **SWISS:** Q99990

Target: VGLL1

Immunogen: KLH conjugated synthetic peptide derived from human VGLL1:

31-130/258.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: VGLL1 is a specific coactivator for the mammalian TEFs. The

mammalian TEF and the Drosophila scalloped genes belong to a conserved family of transcriptional factors that possesses a

TEA/ATTS DNA binding domain.

Applications: WB (1:500-2000)

IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)

Reactivity: (predicted: Human)

Predicted MW.: 29 kDa

Subcellular Nucleus

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.