

bs-9182R

[Primary Antibody]

VGLL1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 51442</p> <p>Target: VGLL1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human VGLL1: 31-130/258.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: VGLL1 is a specific coactivator for the mammalian TEFs. The mammalian TEF and the Drosophila scalloped genes belong to a conserved family of transcriptional factors that possesses a TEA/ATTS DNA binding domain.</p>	<p>Isotype: IgG</p> <p>SWISS: Q99990</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:50-200) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 29 kDa</p> <p>Subcellular Location: Nucleus</p>
---	--	--