

bs-19203R

[Primary Antibody]

NPTII/Neomycin Phosphotransferase 2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

DATASHEET

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: NPTII/Neomycin Phosphotransferase 2

Immunogen: KLH conjugated synthetic peptide derived from Klebsiella pneumoniae Neomycin Phosphotransferase 2: 101-200/264.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: **WB** (1:500-2000)

IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Klebsiella pneumoniae)

Predicted MW.: 32 kDa

SELECTED CITATIONS

- **[IF=2.7]** Wang Jie. et al. A novel truncated variant in SPAST results in spastin accumulation and defects in microtubule dynamics. BMC MED GENOMICS. 2023 Dec;16(1):1-9 WB ;Human. 38066582