

bs-8342R

[Primary Antibody]

CPXM2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 119587	SWISS: Q8N436	Reactivity: (predicted: Human, Mouse, Rat, Pig, Cow)
Target: CPXM2		
Immunogen: KLH conjugated synthetic peptide derived from human CPXM2: 51-150/756.		
Purification: affinity purified by Protein A		Predicted MW.: 83 kDa
Concentration: 1mg/ml		Subcellular Location: Secreted
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: CPXM2 belongs to the peptidase M14 family it contains one F5/8 type C domain and may be involved in cell-cell interactions.		