

bs-8342R

[Primary Antibody]

CPXM2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>Target: CPXM2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human CPXM2: 51-150/756.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: CPXM2 belongs to the peptidase M14 family it contains one F5/8 type C domain and may be involved in cell-cell interactions.</p>	<p>Isotype: IgG</p> <p>GeneID: 119587</p> <p>Applications: WB (1:500-2000) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Cow)</p> <p>Predicted MW: 83 kDa</p> <p>Subcellular Location: Secreted</p>
--	---