

bs-4629R

[Primary Antibody]

Layilin Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 143903</p> <p>Target: Layilin</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human Layilin: 281-382/382.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: LAYN is integral membrane protein, and a hyaluronate receptor. It has been shown to bind specifically to hyaluronan but not to other glycosaminoglycans. It is a recently characterized transmembrane protein with homology to C type lectins, and is a membrane binding site for talin in peripheral ruffles of spreading cells. LAYN may have a role in cell adhesion and motility.</p>	<p>Isotype: IgG</p> <p>SWISS: Q6UX15</p> <p>Applications: WB (1:500-2000) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Cow, Chicken, Dog, Horse)</p> <p>Predicted MW.: 41 kDa</p> <p>Subcellular Location: Cell membrane</p>
--	--