

**bs-17467R**

**[ Primary Antibody ]**

## SHKBP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 92799</p> <p><b>Target:</b> SHKBP1</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human SHKBP1: 31-130/707.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>WB</b> (1:500-2000) <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> Human (predicted: Mouse, Rat, Sheep, Cow, Dog, Horse)</p> <p><b>Predicted MW.:</b> 76 kDa</p> <p><b>Subcellular Location:</b> Cytoplasm</p>
--	---