

---

## Cardiac Troponin I (cTnI) Mouse mAb, AP conjugated

Cat. No:	V3202-AP
Name:	Cardiac Troponin I (cTnI) Mouse mAb, AP conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	17E11K
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Concentration:	>1mg/ml
Buffer:	0.01M PBS(pH7.4) with 0.1% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human Cardiac Troponin I (cTnI)