

**bsm-2024M**

**[ Primary Antibody ]**



www.bioss.com.cn  
sales@bioss.com.cn  
techsupport@bioss.com.cn  
400-901-9800

## HBsAg Mouse mAb

### — DATASHEET —

<b>Host:</b> Mouse <b>Clonality:</b> Monoclonal <b>Target:</b> HBsAg <b>Purification:</b> affinity purified <b>Concentration:</b> 1mg/ml <b>Storage:</b> Size : 50ul/100ul/200ul 0.01M PBS (pH7.4) with 0.02% Proclin300. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. <b>Background:</b>	<b>Isotype:</b> IgM <b>CloneNo.:</b> H2F4	<b>Applications:</b> ELISA (1:5000-10000) <b>Reactivity:</b> (predicted: HBsAg)  <b>Predicted MW.:</b> 44 kDa <b>Subcellular Location:</b> Cell membrane ,Cytoplasm
--	--	---

### — SELECTED CITATIONS —

- **[IF=11.15]** Liu, Chang. et al. Occult infection with hepatitis B virus PreS variants synergistically promotes hepatocellular carcinoma development in a high-fat diet context by generating abnormal ceramides. BMC MED. 2022 Dec;20(1):1-17 IHC ;Mouse. 36058909