

bsm-2024M

[Primary Antibody]



www.bioss.com.cn
sales@bioss.com.cn
techsupport@bioss.com.cn
400-901-9800

HBsAg Mouse mAb

— DATASHEET —

Host: Mouse Clonality: Monoclonal Target: HBsAg Purification: affinity purified Concentration: 1mg/ml Storage: Size : 50ul/100ul/200ul 0.01M PBS (pH7.4) with 0.02% Proclin300. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background:	Isotype: IgM CloneNo.: H2F4	Applications: ELISA (1:5000-10000) Reactivity: (predicted: HBsAg) Predicted MW.: 44 kDa Subcellular Location: Cell membrane ,Cytoplasm
--	--	---

— SELECTED CITATIONS —

- **[IF=11.15]** Liu, Chang. et al. Occult infection with hepatitis B virus PreS variants synergistically promotes hepatocellular carcinoma development in a high-fat diet context by generating abnormal ceramides. BMC MED. 2022 Dec;20(1):1-17 IHC ;Mouse. 36058909