

bsm-2023M

[Primary Antibody]

HBeAg(2E6) Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Monoclonal	CloneNo.: 2E6	Reactivity: (predicted: HBeAg)
Target: HBeAg(2E6)		
Purification: affinity purified by Protein G		
Concentration: 1mg/ml		Predicted MW.: 18 kDa
Storage: Size : 50ul/100ul/200ul 0.01M PBS (pH7.4) with 0.02% Proclin300. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Subcellular Location: Cell membrane ,Cytoplasm
Background:		