bs-16607R

[Primary Antibody]

DILP5 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: DILP5

 $\begin{tabular}{ll} \textbf{Immunogen:} KLH conjugated synthetic peptide derived from Fire fly DILP5: $21-100/108. \end{tabular}$

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) **IF** (1:100-500) ICC/IF (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: Firefly, Fruit Fly)

Predicted 9.5 kDa

Subcellular Location: Secreted