

bs-15515R

[Primary Antibody]

IFIT3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat) Predicted MW.: 56 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 3437	SWISS: O14879	
Target: IFIT3		
Immunogen: KLH conjugated synthetic peptide derived from human IFIT3: 101-200/490.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=4.33]** Tao Yang. et al. Berberine inhibits IFN- γ signaling pathway in DSS-induced ulcerative colitis. SAUDI PHARM J. Saudi Pharm J. 2022 Apr;; WB ;Mouse. 10.1016/j.jsps.2022.03.015
- **[IF=3.448]** Yang T et al. Mechanism of berberine in treating Helicobacter pylori induced chronic atrophic gastritis through IRF8-IFN- γ signaling axis suppressing. Life Sci. 2020 Feb 22;248:117456. WB ;Rat. 32097666