

---

## CYFIP2 Rabbit pAb

Catalog Number: bs-7606R

Target Protein: CYFIP2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Chicken, Dog, Horse)

Predicted MW: 148 kDa

Entrez Gene: 26999

Swiss Prot: Q96F07

Source: KLH conjugated synthetic peptide derived from human CYFIP2: 251-350/1273.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: CYFIP2 is involved in T cell adhesion and p53 dependent induction of apoptosis. It does not bind RNA but is up regulated significantly in CD4+ T lymphocytes from patients with multiple sclerosis. There are 2 isoforms produced by alternative splicing.