
Human VEGF165 Mouse mAb, Biotin conjugated

Cat. No:	V2604-Bio
Name:	Human VEGF165 Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D7
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	Recombinant human VEGF165 Protein