
Anti-Human IgE Mouse mAb, HRP conjugated

Cat. No:	V5402-HRP
Name:	Anti-Human IgE Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5E9
Host/Source:	Mouse
Subtype:	IgG1, Kappa
Reactivity:	Human
Format:	Liquid
Buffer:	lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added.
Applications:	ELISA, ICA
Antigen:	Human IgE