

bs-14764R

[Primary Antibody]

FAM177A1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500)
Clonality: Polyclonal		IHC-F (1:100-500)
GeneID: 283635	SWISS: Q8N128	IF (1:100-500)
Target: FAM177A1		ICC/IF (1:100-500)
Immunogen: KLH conjugated synthetic peptide derived from human FAM177A1: 31-130/213.		Reactivity: (predicted: Human, Mouse, Rat, Pig, Dog)
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		Predicted MW.: 24 kDa
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Subcellular Location: Extracellular matrix, Nucleus

— SELECTED CITATIONS —

- **[IF=5.016]** Huizi Chen. et al. MicroRNA - 7 negatively regulates Toll - like receptor 4 signaling pathway through FAM177A. Immunology. 2021 Jan;162(1):44-57 IF ;Mouse. 32852789