
Cardiac Troponin I (cTnI) Mouse mAb, FITC conjugated

Cat. No:	V3201-FITC
Name:	Cardiac Troponin I (cTnI) Mouse mAb, FITC conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	17C8L
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA Excitation spectrum: 495nm Emission spectrum: 519nm
Antigen:	Recombinant human Cardiac Troponin I (cTnI)

PRODUCT SPECIFIC PUBLICATIONS

[IF=7.923] Panpan Niu. et al. High-sensitive and disposable myocardial infarction biomarker immunosensor with optofluidic microtubule lasing. NANOPHOTONICS-BERLIN. 2022 Jun;; Other ; . 10.1515/nanoph-2022-0208