

bs-0365R-APC

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

human IgG Fc, APC conjugated

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>Target: human IgG Fc</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M PBS(pH7.4). Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.</p>	<p>Isotype: IgG</p>	<p>Applications: IF ICC/IF Excitation Spectrum: 650nm Emission spectrum: 660nm</p> <p>Reactivity: Human</p> <p>Subcellular Location: Secreted</p>
--	----------------------------	--