

human IgG Fc, APC conjugated

Catalog Number: bs-0365R-APC

Target Protein: human IgG Fc

Concentration: 1mg/ml

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IF, ICC/IF

Excitation spectrum: 650nm

Emission spectrum: 660nm

Not yet tested in other applications.

Optimal working dilutions must be determined by the end user.


Reactivity: Human

Storage: 0.01M PBS(pH7.4).


Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.


VALIDATION IMAGES




Sample: Lane 1: Human Serum Primary: Anti-human IgG Fc (bs-0365R) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: kDa Observed band size: 50 kDa



Paraformaldehyde-fixed, paraffin embedded (human tonsil); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (human IgG Fc) Polyclonal Antibody, Unconjugated (bs-0365R) at 1:200 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.



Paraformaldehyde-fixed, paraffin embedded (human colon carcinoma); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (human IgG Fc) Polyclonal Antibody, Unconjugated (bs-0365R) at 1:200 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.



Paraformaldehyde-fixed, paraffin embedded (human colon); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (human IgG Fc) Polyclonal Antibody, Unconjugated (bs-0365R) at 1:200 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.