

bs-0365R-APC

[Secondary Antibodies]

human IgG Fc, APC conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: human IgG Fc Concentration: 1mg/ml Storage: 0.01M PBS(pH7.4). Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.	Isotype: IgG Reactivity: Human Subcellular Location: Secreted
--	--