

---

## D-Dimer Mouse mAb, Biotin conjugated

Cat. No:	V1101-Bio
Name:	D-Dimer Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	G1E9
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	10mM PBS, pH7.4.
Applications:	ELISA
Antigen:	Native D-dimer protein derived from human plasma