
Troponin C (TnC) Mouse mAb, FITC conjugated

Cat. No:	V3401-FITC
Name:	Troponin C (TnC) Mouse mAb, FITC conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	14F11K
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	>1mg/ml
Buffer:	0.01M PBS(pH7.4) with 0.1% Proclin300.
Antigen:	Recombinant human Troponin C (cTnC)