

bs-3134R**[Primary Antibody]****phospho-Ezrin (Thr567) Rabbit pAb****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

— DATASHEET —


Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) Reactivity: Mouse, Rat (predicted: Human, Rabbit, Pig, Cow, Chicken, Dog, GuineaPig, Horse) Predicted MW.: 64 kDa Subcellular Location: Cell membrane ,Cytoplasm
Clonality: Polyclonal		
GeneID: 7430	SWISS: P15311	
Target: Ezrin (Thr567)		
Immunogen: KLH conjugated Synthesised phosphopeptide derived from human Ezrin around the phosphorylation site of Thr567: YK(p-T)LR.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene.		

— VALIDATION IMAGES —

Sample: Small intestine (Mouse) Lysate at 40 ug
Primary: Anti-Phospho-Ezrin (Thr567) (bs-3134R)
at 1/500 dilution Secondary: IRDye800CW Goat
Anti-Rabbit IgG at 1/20000 dilution Predicted
band size: 64 kD Observed band size: 64 kD



Paraformaldehyde-fixed, paraffin embedded
(mouse kidney); Antigen retrieval by boiling in
sodium citrate buffer (pH6.0) for 15min; Block
endogenous peroxidase by 3% hydrogen
peroxide for 20 minutes; Blocking buffer (normal
goat serum) at 37°C for 30min; Antibody
incubation with (Phospho-Ezrin (Thr567))
Polyclonal Antibody, Unconjugated (bs-3134R)
at 1:200 overnight at 4°C, followed by operating
according to SP Kit(Rabbit) (sp-0023)
instructionsand DAB staining.



Paraformaldehyde-fixed, paraffin embedded (rat
kidney); Antigen retrieval by boiling in sodium
citrate buffer (pH6.0) for 15min; Block
endogenous peroxidase by 3% hydrogen
peroxide for 20 minutes; Blocking buffer (normal
goat serum) at 37°C for 30min; Antibody
incubation with (Phospho-Ezrin (Thr567))
Polyclonal Antibody, Unconjugated (bs-3134R)
at 1:200 overnight at 4°C, followed by operating
according to SP Kit(Rabbit) (sp-0023)
instructionsand DAB staining.