

bs-2045R

[Primary Antibody]

alpha-L-arabinofuranosidase Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: alpha-L-arabinofuranosidase

Immunogen: KLH conjugated synthetic peptide derived from Prunus persica alpha-L-arabinofuranosidase: 1-100/677.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: no data available.

Applications: **WB** (1:500-2000)

IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Prunus persica)

Predicted MW.: 74 kDa