

bs-4942R

[Primary Antibody]

HIV1 p55+p24+p17 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: HIV1 p55+p24+p17 Immunogen: KLH conjugated synthetic peptide derived from HIV1 p55+p24+p17: 231-330/497. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	Isotype: IgG Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: HIV) Predicted MW.: 17/24/55/162 kDa
---	---

— SELECTED CITATIONS —

- **[IF=4.307]** Senait Ashenafi. et al. Immunosuppressive Features of the Microenvironment in Lymph Nodes Granulomas from Tuberculosis and HIV-Co-Infected Patients. Am J Pathol. 2022 Jan;; IHC ;Human. 35092727