

bs-19212R

[Primary Antibody]

NEURL1B Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 54492</p> <p>Target: NEURL1B</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human NEURL1B: 1-100/555.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: E3 ubiquitin-protein ligase involved in regulation of the Notch pathway through influencing the stability and activity of several Notch ligands.</p>	<p>Isotype: IgG</p> <p>SWISS: A8MQ27</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Gorilla, Orangutan)</p> <p>Predicted MW.: 59 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
--	--	---