

**bs-12163R**

**[ Primary Antibody ]**

## Calcyphosine 2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Calcyphosine 2

**Immunogen:** KLH conjugated synthetic peptide derived from human Calcyphosine 2/CAPS2: 221-320/557.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

**Reactivity:** (predicted: Human,  
Mouse, Rat, Pig, Sheep,  
Cow, Chicken)

**Predicted  
MW.:** 64 kDa