
beta-Amyloid(1-42) Mouse mAb, Gold conjugated

Cat. No:	V5304-Gold
Name:	beta-Amyloid(1-42) Mouse mAb, Gold conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3D12
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA, ICA
Antigen:	KLH conjugated synthetic peptide derived from human beta-Amyloid(1-42)