

bs-4620R

[Primary Antibody]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

H1N1+H5N1+H3N2 nucleocapsid protein Rabbit pAb

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: H1N1, H5N1, H3N2)
Target: H1N1+H5N1+H3N2 nucleocapsid protein		Predicted MW.: 56 kDa
Immunogen: KLH conjugated synthetic peptide derived from H1N1/H5N1/H3N2 nucleocapsid protein: 151-250/498.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		