

**bs-4621R**

**[ Primary Antibody ]**



www.bioss.com.cn  
sales@bioss.com.cn  
techsupport@bioss.com.cn  
400-901-9800

## **H1N1+H5N1+H3N2 nucleocapsid protein Rabbit pAb**

### **— DATASHEET —**

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>Target:</b> H1N1+H5N1+H3N2 nucleocapsid protein</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from H1N1/H5N1/H3N2 nucleocapsid protein: 201-300/498.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: H1N1, H5N1, H3N2)</p> <p><b>Predicted MW.:</b> 56 kDa</p>
---	--