

bs-6950R

[Primary Antibody]

Hepatitis A virus polyprotein VP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

GeneID: HHAV

Target: Hepatitis A virus polyprotein VP1

Immunogen: KLH conjugated synthetic peptide derived from human Hepatitis A virus polyprotein VP1: 601-700/2227.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Hepatitis A virus (HAV) is classified with the enterovirus group of the Picornaviridae family. Many other picornaviruses cause human disease, including polioviruses, coxsackieviruses, echoviruses, and rhinoviruses (cold viruses). The term hepatitis A (HA) or type A viral hepatitis has replaced all previous designations: infectious hepatitis, epidemic hepatitis, epidemic jaundice, catarrhal jaundice, infectious icterus, Botkins disease, and MS-1 hepatitis.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Human hepatitis A virus)

Predicted MW.: 23/38/245 kDa