

## NSE Mouse mAb

Cat. No:	V1802
Name:	NSE Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6F7
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant human NSE
Pairing Information:	

ELISA		Detection	
		V1801	V1802
Capture	V1801	-	-
	V1802	+++	-
ICA		Detection	
		V1801	V1802
Capture	V1801	-	+++
	V1802	-	-

## PRODUCT SPECIFIC PUBLICATIONS

**[IF=34.9]** Yi Yakun. et al. Nanopore-based enzyme-linked immunosorbent assay for cancer biomarker detection. NAT NANOTECHNOL. 2025 May;;1-8 ELISA ; Human . 40369343