

bs-20851R

[Primary Antibody]

IQCF5 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: IQCF5 Immunogen: KLH conjugated synthetic peptide derived from mouse IQCF5: 81-148/148. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: IQCF5 (IQ Motif Containing F5) is a Protein Coding gene. An important paralog of this gene is IQCF1.	Isotype: IgG Applications: ELISA (1:5000-10000) Reactivity: (predicted: Mouse, Rat) Predicted MW.: 18 kDa Subcellular Location: Cytoplasm
---	--