

NGAL/Lipocalin 2 Mouse mAb

Cat. No: V0605
Name: NGAL/Lipocalin 2 Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 1G7
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant human NGAL protein ([bs-43518P](#))
Pairing Information:

ELISA		Detection						
		V0601	V0602	V0603	V0604	V0605	V0606	V0607
Capture	V0601	-	-	-	-	-	-	-
	V0602	-	-	-	-	-	-	-
	V0603	-	-	-	-	-	-	-
	V0604	-	+	-	-	-	-	-
	V0605	+	+	+++	-	-	-	-
	V0606	+	+	-	-	-	-	-
	V0607	-	-	-	-	-	-	-
ICA		Detection						
		V0601	V0602	V0603	V0604	V0605	V0606	V0607
Capture	V0601	-	-	-	-	-	-	-
	V0602	-	-	-	-	-	-	-
	V0603	-	-	-	-	-	-	-
	V0604	-	-	-	-	-	-	-
	V0605	-	-	-	-	-	-	-
	V0606	+++	+	+	-	-	-	-
	V0607	+	+++	+	-	-	-	-