

MMP9 Mouse mAb

Cat. No: V5901
Name: MMP9 Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 3H9
Host/Source: Mouse
Subtype: IgG2b
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant human MMP9 protein ([bs-43519P](#))
Pairing Information:

ELISA		Detection			
		V5901	V5902	V5903	V5904
Capture	V5901	-	+++	-	-
	V5902	+++	-	-	-
	V5903	+	+	-	-
	V5904	+	+	-	-
ICA		Detection			
		V5901	V5902	V5903	V5904
Capture	V5901	-	-	+++	-
	V5902	-	-	-	-
	V5903	-	-	-	-
	V5904	+	-	-	-