

---

## **AD7c-NTP Mouse mAb, HRP conjugated**

|                |  |
|----------------|--|
| Cat. No:       | V3005-HRP  |
| Name:          | AD7c-NTP Mouse mAb, HRP conjugated   |
| Storage:       | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.     |
| Clone type:    | Monoclonal   |
| Clone No:      | 3B6K   |
| Host/Source:   | Mouse  |
| Subtype:       | IgG1   |
| Reactivity:    | (predicted:&nbsp;Human)  |
| Format:        | Liquid   |
| Concentration: | 1mg/ml   |
| Buffer:        | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.                    |
| Applications:  | ELISA  |
| Antigen:       | KLH conjugated synthetic peptide derived from human Neuronal thread protein AD7c-NTP |