
AD7c-NTP Mouse mAb, HRP conjugated

Cat. No:	V3006-HRP
Name:	AD7c-NTP Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2B4
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human AD7c-NTP Protein