## bs-8547R

## [ Primary Antibody ]

## Bioss

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn

400-901-9800

## **HPV16 L2 Rabbit pAb**

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: HPV16 L2

Immunogen: KLH conjugated synthetic peptide derived from HPV16 L2:

1-100/473.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)

Reactivity: (predicted: HPV16)

Predicted MW.: 52 kDa