

**bs-3202R****[ Primary Antibody ]****phospho-IRS1 (Ser789) Rabbit pAb****BioSS**  
**ANTIBODIES**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**— DATASHEET —****Host:** Rabbit**Isotype:** IgG**Clonality:** Polyclonal**GeneID:** 3667**SWISS:** P35568**Target:** phospho-IRS1 (Ser789)**Immunogen:** KLH conjugated synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Ser789: SS(p-S)SG.**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


**Background:** Insulin receptor substrates (IRS) are responsible for several insulin related activities, such as glucose homeostasis, cell growth, cell transformation, apoptosis and insulin signal transduction. Serine/threonine phosphorylation of IRS1 has been demonstrated to be a negative regulator of insulin signaling and is responsible for its degradation, although IRS1 degradation pathways are not well understood. IRS1 has also been shown to be constitutively activated in cancers such as breast cancer, Wilm's tumors, and adrenal cortical carcinomas, thus making IRS1 phosphorylation and subsequent degradation an attractive therapeutic target. To date there have been four subtypes identified: IRS1, 2, 3 and 4, with IRS1 being widely expressed.**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)


IF (1:100-500)

Flow-Cyt (1µg /test)



ICC/IF (1:100)

**Reactivity:** Human, Rat  
(predicted: Mouse)**Predicted  
MW.:** 132 kDa**Subcellular  
Location:** Cell membrane ,Cytoplasm  
Nucleus**— VALIDATION IMAGES —**

Paraformaldehyde-fixed, paraffin embedded (Rat stomach); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (p-IRS1 (Ser789)) Polyclonal Antibody, Unconjugated (bs-3202R) at 1:400 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.



A549 cell; 4% Paraformaldehyde-fixed; Triton X-100 at room temperature for 20 min; Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min; Antibody incubation with (phospho-IRS1 (Ser789)) polyclonal Antibody, Unconjugated (bs-3202R) 1:100, 90 minutes at 37°C; followed by a conjugated Goat Anti-Rabbit IgG antibody at 37°C for 90 minutes, DAPI (blue, C02-04002) was used to stain the cell nuclei.

Blank control (blue line): HeLa (fixed with 70% ethanol (Overnight at 4°C) and then permeabilized with 0.1% PBS-Tween for 20 min at room temperature). Primary Antibody (green line): Rabbit Anti-phospho-IRS1 (Ser789) antibody (bs-3202R), Dilution: 1µg /10<sup>6</sup> cells; Isotype Control Antibody (orange line): Rabbit IgG . Secondary Antibody (white blue line): Goat anti-rabbit IgG-PE, Dilution: 1µg /test.

Blank control: HepG2. Primary Antibody (green line): Rabbit Anti-phospho-IRS1 (Ser789)

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

antibody (bs-3202R) Dilution: 1µg /10<sup>6</sup> cells;  
Isotype Control Antibody (orange line): Rabbit  
IgG . Secondary Antibody : Goat anti-rabbit IgG-  
AF647 Dilution: 1µg /test. Protocol The cells were  
fixed with 4% PFA (10min at room  
temperature)and then permeabilized with 90%  
ice-cold methanol for 20 min at -20°C. The cells  
were then incubated in 5%BSA to block non-  
specific protein-protein interactions for 30 min  
at room temperature .Cells stained with Primary  
Antibody for 30 min at room temperature. The  
secondary antibody used for 40 min at room  
temperature. Acquisition of 20,000 events was  
performed.