

Recombinant Measles virus hemagglutinin protein, N-His

Catalog Number: bs-42169P

Concentration: >0.5 mg/ml

Species: Measles virus

AA Seq: 59-617/617

Predicted MW: 53.5

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC


Form: Lyophilized or Liquid

Storage: 20mM Tris (pH8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: 麻疹病毒血凝素蛋白(Hemagglutinin protein,H)是介导麻疹病毒吸附、侵入宿主细胞的关键分子,麻疹病毒H蛋白在细胞膜上的表达,以及被麻疹病毒感染的细胞和未感染细胞的直接接触能够引起其受体CD46的表达下调.

VALIDATION IMAGES



The purity of the protein is greater than 80% as determined by reducing SDS-PAGE.