

---

## Recombinant mouse CHAC1 protein, N-His

Catalog Number: bs-42195P

Concentration: >0.5 mg/ml

Species: Mouse

AA Seq: 1-196/223

Predicted MW: 23.8 kDa

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Cation transport regulator like protein 1; Cation transport regulator-like protein 1; ChaC; chac1; CHAC1\_HUMAN.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.