

## Recombinant mouse CD86 protein, C-His (HEK293)

Catalog Number: bs-43590P

Concentration: >0.5mg/ml

Species: Mouse

AA Seq: 24-245/309

Predicted MW: 26.3 kDa

Tags: C-His

Endotoxin: Not tested

Purity: >90% as determined by SDS-PAGE

Purification: AC


Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** This gene encodes a type I membrane protein that is a member of the immunoglobulin superfamily. This protein is expressed by antigen-presenting cells, and it is the ligand for two proteins at the cell surface of T cells, CD28 antigen and cytotoxic T-lymphocyte-associated protein 4. Binding of this protein with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of this protein with cytotoxic T-lymphocyte-associated protein 4 negatively regulates T-cell activation and diminishes the immune response. Alternative splicing results in several transcript variants encoding different isoforms.[provided by RefSeq, May 2011].

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.