

## Human IgE Mouse mAb

Cat. No:	V5404
Name:	Human IgE Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3D3
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted:&nbsp;Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Purified native protein ( <a href="#">bs-2226P</a> )
Pairing Information:	

ELISA		Detection		
		V5402	V5404	V5406
Capture	V5402	-	-	-
	V5404	+++	-	-
	V5406	+++	-	-
ICA		Detection		
		V5402	V5404	V5406
Capture	V5402	-	-	-
	V5404	+++	-	-
	V5406	+++	-	-