

---

## Human IgE Mouse mAb

Cat. No:	V5403
Name:	Human IgE Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1H8
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted:&nbsp;Human)
Buffer:	PBS (pH7.4).
Applications:	ELISA
Antigen:	Human IgE ( <a href="#">bs-2226P</a> )